


VetScan
VS2



THE NEXT GENERATION POINT-OF-CARE ANALYZER



THE NEXT GENERATION POINT-OF-CARE ANALYZER

VetScan VS2 is an innovative state-of-the-art chemistry, electrolytes, immunoassay and blood gas analyzer with uncompromising accuracy from two drops of blood. With VetScan's proven performance, the next generation VetScan VS2 is easier to use than ever before, with an advanced user interface, expanded printing capabilities and precision reference laboratory quality results from 100µl of whole blood. Now you can confidently discuss test results with clients and make confident medical decisions in real time.

FEATURES

- Simple, intuitive, and easy to navigate user interface
- Speed: Results in as few as 10 minutes
- Small platform, space-saving, and elegant design
- Full color touch screen
- Set/Archive reference ranges
- Ability to measure hemolysis, lipemia or icterus
- Software upgrades can be carried out easily via CD drive
- Direct compatibility with supported USB peripherals such as keyboards and printers.
- Stores 5,000 patient records
- VS2 is always ready
- Full range of instrument settings
 - + Five different language choices
 - + Reference range labels
 - + Rotor run status progress bar
 - + Screen saver time options
 - + Back-lit setting options
 - + Sound volume options
 - + Express mode

BENEFITS

- VetScan VS2 system is user-friendly, considerably more cost effective, highly accurate, smaller, lighter, faster and easier to use than ever before.
- Chemistry, electrolyte, immunoassay and blood gas with one simple instrument
- No sample preparation
- Intelligent Quality Control (iQC™)
- Sample flexibility enables testing on multiple species, regardless of size or blood volume, which is vital in smaller patients.
- Complete portability
- Immediate results enables immediate diagnosis and treatment
- Maintenance free
- Every sample is measured for hemolysis, lipemia and icterus
- Complete World-class service, support and warranty



insert sample



insert rotor



read results

| Profile | ALB | ALP | ALT | AMY | AST | BA | BUN | Ca | CHOL | CK | CL- | CRE | GGT | GLOB* | GLU | K ⁺ | Mg | Na ⁺ | PHOS | T ₄ | TBIL | tCO ₂ | TP | UA |
|---|-----|-----|-----|-----|-----|----|-----|----|------|----|-----|-----|-----|-------|-----|----------------|----|-----------------|------|----------------|------|------------------|----|----|
| COMPREHENSIVE DIAGNOSTIC PROFILE #500-0038 | • | • | • | • | | | • | • | | | | • | | • | • | • | | • | • | | • | | • | |
| PREP PROFILE II #500-0026 | | • | • | | | | • | | | | | • | | | • | | | | | | | | • | |
| CRITICAL CARE PLUS #500-0042 | | | • | | | | • | | | | • | • | | | • | • | | • | | | | • | | |
| EQUINE PROFILE PLUS #500-0043 | • | | | | • | | • | • | | • | | • | • | • | • | • | | • | | | • | • | • | |
| LARGE ANIMAL PROFILE #500-0023 | • | • | | | • | | • | • | | • | | • | • | | | | • | | • | | | | • | |
| T ₄ -CHOLESTEROL PROFILE #500-0037 | | | | | | | | | • | | | | | | | | | | | • | | | | |
| AVIAN/REPTILIAN PROFILE PLUS #500-0041 | • | | | | • | • | | • | | • | | | | • | • | • | | • | • | | | | • | • |
| MAMMALIAN LIVER PROFILE #500-0040 | • | • | • | | | • | • | | • | | | | • | | | | | | | | • | | | |

*calculated value

100% rotor test performance guarantee



Abaxis North America
3240 Whipple Road
Union City, California 94587
United States of America
Phone 800 822.2947
Fax 510 441.6150

Abaxis Europe
Heidelberger Landstrasse 230
D-64297 Darmstadt
Germany
Phone +49 6151 350 79 0
Fax +49 6151 350 79 11