

VetScan<sup>5</sup>  
HM5



STATE-OF-THE-ART HEMATOLOGY  
5-PART DIFFERENTIAL AT THE TOUCH OF A BUTTON

# VetScan<sup>5</sup> HM5

## STATE-OF-THE-ART HEMATOLOGY 5-PART DIFFERENTIAL AT THE TOUCH OF A BUTTON

VetScan HM5, a fully-automated 5-part cell counter specifically designed for veterinary applications, offers a comprehensive 22-parameter complete blood count with cellular histograms in just minutes. Its superior performance, elegant design, ease of use, true database management capability, and minimal maintenance make it the optimal hematology system for veterinary clinics, research laboratories, and pharmaceutical/biotech companies.

### FEATURES

- 50 µl sample size
- 5-part differential on Cat, Dog, Horse
- Seven other species available in 3-part differential mode
- Small, space-saving, and elegant design with large LCD display
- Pre-dilute mode available for challenging samples
- Simple, intuitive, and easy to navigate software interface
- 22 parameters, including eosinophil counts and eosinophil percentage
- Advanced and integrated self-cleaning system
- Large database storage capacity, simplicity, flexibility, and usability
- CD drive and USB ports
- Integrates to the VetScan and data management systems
- User-reminder functions
- Multi-user mode

### BENEFITS

- 10 species available now
- Small footprint
- Easy to use
- Reliable results, including out of range values
- Rapid test turn-around with minimal hands-on time
- Customizable reference ranges
- Optimized performance with minimal maintenance
- Automatically saves up to 2000 test records. All downloadable through USB pen drive and integrated data management system
- Allows for software updates and future instrument upgrades
- Integrated result reporting and comprehensive laboratory analysis
- Software tracks reagent usage and notifies user to replace reagents when necessary
- Eco-friendly reagents

### THREE EASY STEPS

Seconds of hands-on time is needed for performing CBC measurements using the VetScan HM5. All you need to do is 1) enter patient information, 2) insert a well-mixed blood sample and press start, 3) read results. The test results and four cellular histograms will be displayed in minutes.




enter information



insert sample and press start



read results



# The all-new VetScan® HM5. It's reproducible.

re-pro-duc-ible (rē'prə-dōōs-sə-bəl) *adj.* it's results are capable of being produced again, it's able to be recreated; it's on the spot with machine-like precision; it's the difference between hitting the dart board and hitting a bull's-eye; it's bowling a 300 or batting a 1000; it's a carbon copy, an exact duplicate, a repeat performance; it's consistently accurate; it's 5-part differential over and over again; most importantly, it's the difference between a confident diagnosis and an educated guess.

The all-new VetScan HM5. You can count on it.



## The all-new VetScan<sup>®</sup> HM5. It's accurate.

ac·cu·rate (ak'yər ət) *adj.* it's careful and exact; it's correct; without mistakes or errors; it's the truth, the whole truth and nothing but the truth; it's reference lab quality; it's the skinny, the story, the fine print, the real deal; it's 5-part differential with unquestionable integrity; it's the difference between cold, hard facts and just plain wrong; it's a true representation; it's on the money, on the spot, right on target; it's the difference between a confident diagnosis and an educated guess. The all-new VetScan HM5. You can count on it.

