

**mindray**

## **DP-50 Vet**

Digital Ultrasonic Diagnostic Imaging System

# The Premium System for Your Choice



# DP-50 Vet

With comprehensive configurations and an integrated design, NEW DP-50 Vet is the outcome of Mindray's continuous and determined efforts towards making veterinary healthcare more efficient, effective and accessible for all. The DP-50 Vet, a premium Black & White ultrasound system, offers more than your expectation. With a smart new shape and full screen, enhanced mobility and more convenient operations, this reliable and affordable system will enhance general practices and is well-suited for multiple veterinary clinical needs.

## Enhanced Veterinary Diagnostic Confidence

- **Power Doppler:** Identify blood flow under various veterinary clinical settings
- **PW Doppler and Auto Trace:** Reveal details of blood flow for more comprehensive diagnosis
- **Veterinary Applications:** Measurement, Comment, Report and Body mark designed for canines, felines, equines, bovines, ovines and more



## Mobile Design for Veterinary Daily Use

- Sleek, streamlined, and compact shape
- 60 degree tilting angle adjustable monitor
- 2 universal transducer connectors
- 500GB hard disk for large patient data storage
- 1.5h scanning with rechargeable battery
- Slim trolley for easy mobility
- Convenient hand-carried bag for easy transportation

## Dedicated Veterinary Workflow

- **iTouch™:** Quick image optimization in one touch
- **iZoom™:** Instant full screen image expansion with one button
- **iStorage™:** Directly transfer images and reports to PC via network cable
- **User-defined keys:** programmable personalized operations tailored for veterinarian



Canine Gallbladder



Canine Bladder

P.I. La Barrena. Canal Imperial de Aragón S/N.  
31500 TUDELA. NAVARRA. (SPAIN)  
Tel.: 948 82 46 55 • 948 82 17 13  
e-mail: [cvm@cvm.es](mailto:cvm@cvm.es) • web: [www.cvm.es](http://www.cvm.es)

